D3.2 SOPs (standard operation procedures) for samples and data flow among surveillance sites and national PH labs

Work Package 3: Enhancement & Consolidation of WGS- & PCR-based Methods for Public Health Action

Author(s): Grigorios Spanakos, Ioanna Spiliopoulou, Anastasia Flountzi

This project has received funding from the European Union's Health and Digital Executive Agency (HaDEA) under Grant Agreement No 101102440.

Disclaimer

Funded by the European Union. Views and opinions expressed are however those of the author(s) only and do not necessarily reflect those of the European Union or HaDEA. Neither the European Union nor the granting authority can be held responsible for them.

Copyright message

© HERA 2 Consortium, 2024

This deliverable contains original unpublished work except where clearly indicated otherwise. Acknowledgement of previously published material and of the work of others has been made through appropriate citation, quotation or both. Reproduction is authorized provided the source is acknowledged.

Document Information

Grant Agreement Number	101102440	Acronym	HERA 2	
Full Title	Consolidation of WGS/RT-PCR and infrastructure processes in surveillance and outbreak investigation activities			
Call	EU4H-2022-DGA-MS-IBA-1			
Topic	EU4H-2022-DGA-MS-IBA-01-02			
Type of Action	EU4H-PJG			
Start Date	01.10.2022	Duration (in months)	42	
HaDEA Project Officer	Marie de Looz-Corswarem			
Task	T3.1			
Work Package	3			
Date of Delivery	Contractual	M15/Dec 2023	Actual	M16/Jan 2024
Nature	R – document, report	Dissemination Level	PU - Public	
Lead Beneficiary	EODY			
Lead Author	Grigorios Spanakos	Organization	EODY	
Other authors	Ioanna Spiliopoulou, Anastasia Flountzi			
Reviewer(s)	Daniel Polzer, Julia Steger			

Document History

Version	Issue Date	Stage	Changes	Contributor
0.1	02.10.2023	First draft generated in Zadar steering meeting		Whole consortium
0.2	19.12.2023	Draft shared with consortium		GS, IS, AF
0.3	16.01.2024	Updated draft with translations into all national languages for proofreading by		GS
1.0	31.01.2024	Revised and quality-checked draft		DP, JS

1. EXECUTIVE SUMMARY	5
2. ABBREVIATIONS AND ACRONYMS	5
3. BACKGROUND INFORMATION	6
4. ENVIRONMENTAL SAMPLES	7
5. HUMAN AND ANIMAL SAMPLES	7
6. METADATA.....	7
7. STATUS OF CONSORTIUM PARTNERS	9
8. DISCUSSION	10
9. SOP – Clinical (and veterinary) samples	11
10. SOP - Water samples.....	15
11. SOP - Food Samples	18
12. German Translation	22
12.1 SOP - Klinische (und tierärztliche) Proben.....	22
12.2 SOP - Wasserproben	26
12.3 SOP - Lebensmittelproben	29
13. Croatian translation	33
13.1 SOP – humani i veterinarski (životinjski) uzorci	33
13.2 SOP – uzorci vode.....	37
13.3 SOP – Uzorci hrane	40
14. Hungarian translation	44
14.1 SOP - Klinikai (és állatgyógyászati) minták.....	44
14.2 SOP - Vízminták.....	48
14.3 SOP - Élelmiszerminták	51
15. Greek translation.....	56
15.1 SOP – Κλινικά (και κτηνιατρικά) δείγματα	56
15.2 SOP – Δείγματα ύδατος.....	60
15.3 SOP – Δείγματα τροφής	63

1. EXECUTIVE SUMMARY

Sample analysis may be divided to preanalytical, analytical and post analytical stages. The preanalytical stage is of major importance as errors during this stage affect the validity of the results, even if subsequent stages were performed correctly. The preanalytical stage includes sampling and transport of the sample and metadata to the laboratory. Sampling and transport are described in detail in various documents published by the ISO, or National and International Health agencies, like World Health Organization (WHO).

Molecular Laboratories of Public Health Institutions in the context of the One Health approach, have to cope with a variety of sample types and aims. Samples may be of environmental, food, animal, and human origin, while the aim could be epidemiologic surveillance, outbreak investigation or diagnosis.

Environmental samples such as water, wastewater, and food, usually require the application of methods that increase the pathogen's concentration, and this fact defines transport conditions. Official procedures are described in International Organization for Standardization (ISO) documents. Clinical and veterinary samples are usually used directly for nucleic acid extraction, and freezing is considered the best method of preservation during transport to the laboratory. In practice however less demanding alternative methods are usually used.

Metadata are information associated with characteristics of the sample. Metadata make the sample and the data derived from the sample useful. In the context of Public Health, metadata may be useful to the laboratory or to the epidemiology department. Clinical specimens' metadata may include personal data, and in certain cases this is also true for veterinary specimens. The existence of personal data poses a problem regarding the transfer of the information from the person collecting the specimen to the laboratory and finally to the epidemiology department, as personal data should not become publicly available.

To fulfil the requirements for General Data Protection Regulation (GDPR) compliance only trained and authorised personnel should have access to personal data. This can be attained using database linked to a web platform for data entry, and data will be accessible only to predefined competent personnel. As all consortium partners do not currently have such a suitable infrastructure, GDPR compliance can be achieved with trained authorized personnel and encryption for data transfer.

2. ABBREVIATIONS AND ACRONYMS

CDC	Centers for Disease Control and Prevention
GDPR	General Data Protection Regulation
ISO	International Organization for Standardization
LIMS	Laboratory Information Management Systems
SOPs	Standard Operating Procedures
WHO	World Health Organization

3. BACKGROUND INFORMATION

Standard Operating Procedures (SOPs) are written step-by-step instructions to perform a task. SOPs aim to avoid ununiform performance and errors during the task execution. SOPs should be clear, comprehensive, and well written, easy to understand and follow and finally applicable and up to date.

Public Health Laboratories, among other tasks, detect, analyse, and characterize pathogenic microorganisms in various types of samples, originating from patients, animals, or the environment, aiming to public health protection from infectious agents that may pose threat to people's health. Microbiological results are important for diagnosis and treatment of patients, surveillance, and intervention. Nowadays, pathogenic microorganisms "travel" across large distances, with the most recent example being the SARS-CoV-2 coronavirus that originated from China and spread all over the world causing, in a short time, the COVID-19 pandemic. It is therefore necessary to ensure that Public Health laboratory results are reliable, accurate, consistent, and reproducible across laboratories.

The activities of a Microbiological laboratory can be divided in three stages:

- the preanalytical step comprising sampling, transportation, and storage until analysis
- the analytical stage that includes the analytical process
- the postanalytical stage that deals with reporting and interpretation

While every stage is important for results accuracy, the preanalytical stage directly affects all subsequent procedures. The preanalytical stage includes two parts; one is the sample collection and handling up to introduction to the analytical stage; the other is the collection, storage, and transfer to the competent authority, of the information that is related to the characteristics of the sample (metadata) which are necessary for the interpretation and evaluation of the results produced during the analytical stage. Wrong metadata or wrong sampling or even wrong type of sample lowers the validity and usefulness of the results produced during the subsequent stages.

Public Health Laboratories performing molecular detection and NGS under the One Health approach have to deal with several different types of samples. One crucial factor for differentiation of the sample types is the sample's origin, that affects the sampling procedure, transport conditions, and required metadata. Depending on the origin, samples may be environmental (water, wastewater, food, feed), animal or human. Additionally, for all types of samples, nucleic acid isolation for NGS applications may be performed directly on the sample or on bacterial colonies derived from culture of the sample and this is affecting the transport and storage conditions required.

An important aspect of the preanalytical stage is the preservation of the sample. All samples are biologically active and removal of the sample from the natural environment results in changes in the microbial composition, as microorganisms could die or grow depending on the conditions. To avoid such changes, the sample should ideally be processed immediately after collection, but this is seldom the case for Public Health laboratories. To maintain the original composition of the sample, specific measures should be taken depending on the pathogen(s) of interest. Chemical preservation, which is common in classical clinical microbiology practices, is not usable for molecular analysis as chemicals render nucleic acids useless.

Besides the type of samples, the analysis's purpose also adds another degree of complexity to the information (metadata) provided with the sample. Simple surveillance population studies do not always require personal information of the participants, while the opposite is true for outbreak investigation, that may require detailed health status reports and communication details. Also, depending on the pathogen, different information is necessary. For example, for STDs information on sexual behaviour is needed, while for COVID-19 it is not.

4. ENVIRONMENTAL SAMPLES

Environmental samples usually are complex and have low numbers of pathogens per unit of mass or volume. To facilitate detection, samples are subjected to various procedures to increase the concentration and/or remove unwanted materials (e.g., centrifugation, filtration, pre-enrichment, etc). The efficiency of these procedures depends on the integrity of the pathogens' structure, consequently such samples should be maintained until analysis in a way that preserves the structure of the microbes. The usual way is to refrigerate the sample and keep it at a low temperature. Such samples should not be frozen, unless they already are e.g., frozen food.

Procedures of environmental sampling are included in various ISO methods, which are official guidelines therefore they should be considered as the best approach for sampling, transfer, storing and processing of these types of samples. ISO 19458:2006 describes sampling of water for microbiological analysis, while ISO/TS 17728:2015 describes sampling of food or feed for microbiological analysis. Furthermore, ISO 18593:2018 specifies horizontal methods for surface sampling.

5. HUMAN AND ANIMAL SAMPLES

Samples such as blood, urine, stool, or tissue, are usually subjected directly to nucleic acid isolation for application of molecular methods, therefore it is recommended to be transported frozen (-20 °C or lower) to the laboratory. Frozen samples should not melt until the extraction of DNA and/or RNA, as repeated cycles of freezing and melting result in breaking of nucleic acid molecules.

Frozen samples, however, require dry ice and special packaging, which are not always available. Alternative specific transport mediums may be used, and refrigeration to inhibit growth, until arrival to the laboratory. Transport media are buffered solutions containing minimum nutrients that support viability of the microbe of interest and does not favour growth of others. If bacterial culture before molecular analysis is needed, samples should not be frozen. They should be refrigerated and maintained at a low temperature.

WHO (WHO. Laboratory testing for coronavirus disease (COVID-19) in suspected human cases. Interim guidance 19 March 2020), CDC and other Health Agencies have published guidelines for collection, transport, and storage of human specimens, such as blood, stool, nasopharyngeal swabs etc.

6. METADATA

Metadata are information relevant to the sample that answer basic questions about the sample, like "when, who, why," and more, that are registered during the whole analysis procedure by all persons dealing with the sample, such as the person who collected the sample, the person that received the sample, the analysts, and the laboratory supervisor. Metadata in the context of Public Health authorities can be divided into those necessary for the laboratory and those necessary for epidemiology. Table 1 presents laboratory useful metadata for different sample types.

Environmental sample's metadata do not include information that could reveal personal data. On the contrary, clinical sample's metadata usually include personal information, while veterinary sample's metadata may include personal data such as the animal owner's name. The inclusion of personal data could pose a problem when communicating the information to all involved stakeholders of the Institution, as this information should not be publicly available due to restrictions imposed by the GDPR. Consequently, specific measures should be applied. Indicative metadata for animal and human samples are presented in Table 2.

As already mentioned, the list of metadata varies depending on the purpose of the analysis and the pathogen.

Table 1. Necessary laboratory metadata until analysis for environmental, animal, and clinical samples.

	Environmental samples	Animal samples	Clinical samples
1	Sample ID or Barcode	Sample ID or Barcode	Sample ID or Barcode
2	Date and time of sampling	Date and time of sampling	Date and time of sampling
3	Type of sample (e.g. wastewater, milk, lake water, meat etc)	Type of sample (e.g. urine, stool, biopsy etc)	Type of sample (e.g. blood, nasopharyngeal swab, etc)
4	Additional sample information	Additional sample information	Additional sample information
5	List of parameters to analyze	List of parameters to analyze	List of parameters to analyze
6	Location of sample collection	Location of sampling or owner or farm address	Patient name and address
7	Date and time of arrival at the laboratory	Date and time of arrival at the laboratory	Date and time of arrival at the laboratory
8	Sample Temperature	Sample Temperature	Sample Temperature
9	Reasons to reject the sample	Reasons to reject the sample	Reasons to reject the sample
10	Laboratory ID	Laboratory ID	Laboratory ID
11	Person that received the sample	Person that received the sample	Person that received the sample
12	Date and time that laboratory analysis started	Date and time that laboratory analysis started	Date and time that laboratory analysis started

Table 2. Animal and human samples' indicative metadata.

Veterinarian sample	Human sample
Animal species	Patient name
Animal age or sex	Patient age and sex
Animal symptoms (if any)	Patient symptoms (if any)
Onset of symptoms (if applicable)	Onset of symptoms (if applicable)
Travel history	Travel history
Vaccination status (if applicable) (for serological tests)	Vaccination status (if applicable)
	Co-morbidities
Medication (if applicable)	Medication (if applicable)
Additional information on the animal	Additional information on the patient

7. STATUS OF CONSORTIUM PARTNERS

Partner's laboratories responses regarding SOPs and accreditations were as follows:

The department of molecular biology of the veterinary department in AGES Austria is accredited according to ISO 17025 for PCR. There is a SOP regarding sample logistics within the institute in the veterinary department.

The National Public Health Center (NCPHP) of Hungary is accredited by the National Accreditation Authority according to the current ISO 15189 and ISO 17025 standards, regarding NGS, and only *Listeria* WGS is involved in this accreditation period. There is a limited number of SOPs.

The Genotyping Department (Division for Microbiology) of Croatian Institute of Public Health is not accredited and is currently in the process of writing SOPs.

In Greece none of the Public Health laboratories are accredited for WGS. Only the Regional Public Health Laboratory of Thessaly is accredited according to ISO 17025 for detection of SARS-CoV-2 using real time PCR.

Table 3. Metadata used by the consortium partners.

Austria	Hungary	Croatia	Greece
Sample provider (name, address)	Patient name		Patient name
Laboratory ID	Laboratory ID	Laboratory number	Laboratory ID
Type of sample (tissue, blood, swab, semen)	Type of sample	Type of sample (clinical sample/bacterial isolate/DNA or RNA extract)	Type of sample
about the sample species (pig, chicken....) Host	Host	Host	
Collection date	Collection date	Sampling date	Sampling date
Sampling location (in case of wild animals gps coordinates, in case of domestic animals the farm ID)	Location (just country)	Location (county)	Location
	Gender		Gender
	Patient status if possible, travel history, animal contact, symptoms, vaccination status (if applicable)		
Requested analysis			Requested analysis

All consortium partners reported that they use excel files for data storage, while Austria uses LIMS for data storage. Samples are accompanied by a form filled with sample ID and metadata. Laboratory personnel fill in spreadsheet files with samples' metadata along with the laboratory results and send them by mail to the epidemiology department, or through LIMS.

In Greece, for COVID-19 samples' entry, a web platform is used, and sample barcode and metadata are recorded on site, stored in a database and are immediately available to the Directorate of Epidemiology. Laboratory submits the results through the same web platform, using only the barcode of the sample. In that way the laboratory does not have any information that could infringe the GDPR. This approach resulted in lower turnaround times from sample entry at the laboratory to results and lower laboratory workload.

8. DISCUSSION

Ideally the laboratory should only receive metadata required for the analysis, so that analysis is objective and GDPR is not violated. This would be possible using a database linked to a web platform accessible for data entry to all parties involved, i.e. physicians, veterinarians, samplers, epidemiologists, and laboratory personnel. The person collecting the sample enters all required data (including personal data), but all other parties may access only a restricted set of data relevant to their role in the procedure. Any data that could infringe the GDPR should be accessible only to trained and authorized personnel. The laboratory supervisor, for example, who is responsible for the final approval for results release should have access to the clinical condition of the patient, for evaluation of the results, as in certain cases additional analysis or repetition of the test may be necessary. As mentioned above a similar approach was applied successfully in Greece.

Currently all consortium partners do not have suitable databases for remote data entry but use spreadsheets, while Austria uses additionally the Ridom SeqSphere+ Software and LIMS software for local data management. A possible scenario for data communication among involved parties, i.e. person performing sampling, laboratory, and department of epidemiology, with respect to GDPR would be:

The person who collects the sample, fills in the form with the necessary data (either paper or spreadsheet).

Depending on the form, transport along with the sample, or send (encrypted and secured with a strong password) by email to the laboratory. At the laboratory, authorised personnel, may also import the data into LIMS software if available.

Only authorized laboratory personnel process the data, adds the results, and sends the file to the department of epidemiology encrypted and secured as above, or data are disseminated through LIMS software.

9. SOP – Clinical (and veterinary) samples

Clinical and veterinary samples' testing is performed at Public Health Laboratories for surveillance, outbreak investigation, and diagnosis.

Collection of clinical samples is described in various documents such as “J. Michael Miller, Ph.D. Handbook of Specimen Collection and Handling in Microbiology. Second Edition. 1985”, “WHO/CDS/CSR/EDC/2000.4. Guidelines for the collection of clinical specimens during field investigation of outbreaks (available at: https://iris.who.int/bitstream/handle/10665/66348/WHO_CDS_CSR_EDC_2000.4.pdf), “Protocol to investigate non-seasonal influenza and other emerging acute respiratory diseases. Geneva: World Health Organization; 2018. (available at: <https://www.who.int/publications/i/item/WHO-WHE-IHM-GIP-2018.2>)” and “Koneman's Color Atlas and Textbook of Diagnostic Microbiology 6th Edition”.

Veterinary sample collection is described in “Bildfell R. Collection and submission of laboratory samples from animals. Merck Veterinary Manual” (available at: merckvetmanual.com/clinical-pathology-and-procedures/collection-and-submission-of-laboratory-samples/collection-and-submission-of-laboratory-samples-from-animals), accessed November 2023.

Clinical samples should be collected by a physician (or a veterinarian in case of animals), or other competent and authorized personnel. Samples should be placed in suitable containers with size proportional to the sample size. They should be identified with a properly marked strip of masking tape, or a barcode strip.

The person responsible for sampling, should fill in the relevant data fields of the sample form (Annexes I and II for clinical (human) and veterinary samples respectively).

If the sample form is in a paper form, it is sealed in a plastic bag and inserted in the sample packaging. If it is a spreadsheet it is sent, encrypted, and secured by a strong password, by email to the receiving laboratory.

In general, specimens for virus detection should reach the laboratory as soon as possible after collection. For most clinical samples, the recommended shipment time is maximum 2 days at 2-8 °C. If delivery will be delayed for more than the recommended shipment time of the specific specimen type, specimen should be frozen at -70 °C. In case the samples will be directly subjected to DNA and/or RNA extraction without prior processing e.g., centrifugation, they should be sent to the laboratory frozen at -70 °C (dry ice), when possible. In case bacteriological analysis (culture) is required before the application of molecular methods, the sample must be delivered at the laboratory within 24 hours at 4 °C for stool or lower respiratory specimens, or at room temperature for rectal swabs and upper respiratory specimens.

(WHO Laboratory testing for coronavirus disease (COVID-19) in suspected human cases: interim guidance (<https://www.who.int/publications/i/item/10665-331501>).

If delivery time is expected to be less than the maximum acceptable, the sample is placed on refrigerant gel packs, or at 2-8 °C until shipment, unless it is otherwise indicated.

The sample is packed in suitable container that maintains the recommended temperature (2-8 °C or -70 °C) and is in accordance with the local/international regulations about transportation of hazardous material class B.

Upon arrival at the laboratory, an analyst who is authorized to access personal data, records the temperature, and notes its physical condition. Accordingly fills in the relevant fields in the sample form.

Samples are stored at appropriate conditions, as indicated, until analysis.

The sample is rejected if:

Samples are unlabelled or mislabelled.

Sample temperature is not the indicated with respect to the sample type and requested test.

Length of time between sample collection and receipt is longer than the recommended

Sample container leaks

There are visible cracks or any form of damage to the container

Sample volume is less than the recommended

Mandatory fields in the data form are not filled in

If rejected, in the "Reasons to reject the sample" field of the form the exact reason for rejection should be mentioned.

The laboratory analyst (authorized to access personal data) completes the "Date and time that laboratory analysis started" field when starting the analysis procedure. If LIMS software is available all metadata are imported.

When results are ready, the form and the results report are sent by email, encrypted, and secured by a strong password, to the epidemiology department, or data may be accessible through LIMS software.

ANNEX I

Metadata for clinical (human) samples

	Clinical sample
1	Sample ID or Barcode*
2	Date and time of sampling*
3	Type of sample (e.g. blood, nasopharyngeal swab, etc)*
4	Sample shipment condition (frozen, refrigerated, room temperature)*
5	Location of sample collection (country only)*
6	List of parameters to analyze*
7	Patient name and address
8	Date and time of arrival at the laboratory*
9	Sample temperature*
10	Reasons to reject the sample (Yes/No, description)*
11	Laboratory ID*
12	Name of the person that received the sample*
13	Date and time that laboratory analysis started*
14	Location of sample collection (exact)
15	Patient name
16	Patient age
17	Patient sex
18	Patient symptoms (if any)
19	Onset of symptoms (if applicable)
20	Travel history
21	Vaccination status (if applicable)
22	Co-morbidities
23	Medication (if applicable)
24	Additional information on the patient

Data fields 1-7 and 14-24 are completed by the person performing the sampling, while 8-14 are completed by the laboratory personnel (fields marked with asterisk are mandatory).

Note that patient data fields (14-24) are subjected to changes depending on the pathogen and the purpose of the analysis. Also, fields that are not indicated as mandatory, may become mandatory if for example “contact tracing” of the patient is required.

ANNEX II

Metadata for veterinary samples

	Veterinary sample
1	Sample ID or Barcode*
2	Date and time of sampling*
3	Animal species*
4	Type of sample (e.g. urine, stool, biopsy etc)*
5	Additional sample information
6	List of parameters to analyze*
7	Location of sampling or farm address*
8	Date and time of arrival at the laboratory*
9	Sample temperature*
10	Reasons to reject the sample*
11	Laboratory ID*
12	Name of the person that received the sample*
13	Date and time that laboratory analysis started*
14	Animal age or sex
15	Animal symptoms (if any)
16	Onset of symptoms (if applicable)
17	Medication (if applicable)
18	Additional information on the animal

Data fields 1-8 and 14-18 are completed by the person performing the sampling, while 8-13 are completed by the laboratory personnel (fields marked with asterisk are mandatory).

Note that fields (14-24) are subjected to changes depending on the pathogen and the purpose of the analysis.

10. SOP - Water samples

Water from various sources and wastewater samples are subjected to regular testing to assess quality parameters and presence of pathogens for Public Health surveillance. They are also useful for detecting the source of a waterborne outbreak, measuring the pathogens' load in a community for epidemiological reasons, or detecting antibiotic resistance markers in animal husbandry farms, etc.

Sampling should be performed according to ISO 19458:2006 and/or ISO methods corresponding to specific pathogens listed below (Annex I). Wastewater based SARS-CoV-2 surveillance should be performed according to the WHO guidelines available in <https://iris.who.int/bitstream/handle/10665/372995/9789240080638-eng.pdf?sequence=1> (accessed 2023-12-10)

Water samples should be collected in clean sterile bottles, that may be made either from glass, being reusable after sterilization, or from plastic material. Whatever the material, bottles intended for sampling water for human consumption should contain sodium thiosulfate (20mg/L) to neutralize residual chlorine.

If a large volume of water is required (e.g. for *Cryptosporidium*, Norovirus), larger containers may be used, but they should always be clean and sterile. Alternatively, filtration or sedimentation may be performed in situ, to avoid complexities related with the transport and refrigeration of large volumes of liquids.

If sampling should be performed by submerging the bottle, it should be sterile both inside and outside, and kept in a suitable containment, e.g. in a plastic bag or in aluminum foil until used for sampling. For submergence, the bottle should be held with a sterile glove, or a sterile tool (e.g. sterilizable poles or forceps).

Finally, sampling can be performed using water pumps or special samplers suitable for the sampling purpose. Regardless of the sampling method applied, contamination of the sample from the sampling equipment or the environment must be avoided.

After sampling the bottle is sealed and the cap is wrapped with aluminum foil or parafilm to prevent contamination.

The sample should be identified with a properly marked strip of masking tape, or a barcode strip.

The person that performed the sampling fills in the sample form in the relevant data fields (Annex II).

If it is in a paper form, it is sealed in a plastic bag and inserted in the sample packaging. If the form is a spreadsheet it is sent by email to the receiving laboratory.

Samples should be sent to the laboratory refrigerated, and sample temperature ideally should be maintained at 5 ± 3 °C until analysis, unless other conditions are indicated for the specific pathogen.

Upon arrival at the lab, the person who receives the sample records the temperature and notes its physical condition. Accordingly fills in the relevant fields in the sample form. Samples are stored at appropriate conditions, as indicated.

The sample is rejected if:

Samples are unlabelled or mislabelled.

Sample temperature is over 10 °C (unless it is differently specified in the corresponding ISO method)

Time from sampling is longer than the recommended

Sample container leaks

There are visible cracks or any form of damage to the containers

Sample volume is less than the recommended

Mandatory fields in the data form are not filled in

If rejected, in the “Reasons to reject the sample” field of the form the exact reason for rejection should be mentioned.

The laboratory analyst completes the "Date and time that laboratory analysis started" field when starting the analysis procedure. If LIMS software is available all metadata are imported.

When results are ready, the form is sent by email to the epidemiology department along with the results report, or data may be accessible through LIMS software.

ANNEX I

Official guidelines for pathogen detection and/or enumeration.

For bacterial colony isolation to perform WGS:

ISO 11731:2017. Water quality. Enumeration of *Legionella*

ISO 7899-2:2000 Water quality — Detection and enumeration of intestinal enterococci — Part 2: Membrane filtration method

ISO 19250:2010 Water quality — Detection of *Salmonella* spp.

ISO 14189:2013 Water quality — Enumeration of *Clostridium perfringens* — Method using membrane filtration

ISO 9308-3:1998/Cor 1:2000 Water quality — Detection and enumeration of *Escherichia coli* and coliform bacteria — Part 3: Miniaturized method (Most Probable Number) for the detection and enumeration of *E. coli* in surface and wastewater — Technical Corrigendum 1

ISO 9308-1:2014/Amd 1:2016 Water quality — Enumeration of *Escherichia coli* and coliform bacteria — Part 1: Membrane filtration method for waters with low bacterial background flora — Amendment 1

ISO 6461-2:1986 Water quality — Detection and enumeration of the spores of sulfite-reducing anaerobes (clostridia) — Part 2: Method by membrane filtration

For molecular detection of pathogens by qPCR:

ISO/TS 12869:2019 Water quality. Detection and quantification of *Legionella* spp. and/or *Legionella pneumophila* by concentration and genic amplification by quantitative polymerase chain reaction (qPCR).

“Environmental surveillance for SARS-CoV-2 to complement other public health surveillance” The most recent guidelines for SARS-CoV-2 detection and enumeration in Wastewater samples (<https://iris.who.int/bitstream/handle/10665/372995/9789240080638-eng.pdf?sequence=1>). The procedure may be used for surveillance of other viruses shed in feces, such as influenza viruses, RSV.

As official molecular detection methods are limited, concentration methods described in ISO documents may be used such as ISO 15553:2006 (Water quality — Isolation and identification of *Cryptosporidium* oocysts and *Giardia* cysts from water) and nucleic acids are isolated from the concentrate for molecular detection.

For metagenomic analysis, and/or WGS by selective whole genome amplification, DNA and/or RNA may be isolated from sample concentrates, or directly from the sample.

ANNEX II

Metadata for water samples

	Water or wastewater sample
1	Sample ID or Barcode*
2	Date and time of sampling*
3	Type of sample (e.g. wastewater, lake water, tap water, etc)*
4	Additional sample information (surface sampling, sampling depth, etc)
5	List of parameters to analyze*
6	Location of sample collection*
7	Link to outbreak (if it is applicable)
8	Date and time of arrival at the laboratory*
9	Sample temperature*
10	Reasons to reject the sample (Yes/No, description)*
11	Other remarks (if any)
12	Laboratory ID*
13	Name of the person that received the sample*
14	Date and time that laboratory analysis started*

Data fields 1-7 are completed by the person performing the sampling, while 8-14 are completed by the laboratory personnel (fields marked with asterisk are mandatory).

11. SOP - Food Samples

Food samples are tested regularly for Public Health surveillance and/or during foodborne outbreak investigations. Sampling must be performed according to ISO/TS 17728:2015 and/or ISO methods corresponding to specific pathogens listed below (Annex I).

Samples should be properly collected and handled and be representative of the sampled lot. Sampling procedures must be applied uniformly. A representative sample is essential when pathogens or toxins are sporadically distributed within the food. The number of units that comprise a representative sample from a designated lot of a food product must be statistically significant.

Sample units should consist of at least 100gr of the sample whether alone or split in parts, e.g., four 25 g containers could constitute a sample unit.

Sample units should be collected randomly to ensure that a sample is representative of the lot. Empty sample containers should be used as controls when containers are used to collect samples.

Samples should be sent to the laboratory in the original unopened containers, whenever this is possible. If not possible, representative portions should be transferred to sterile containers under aseptic conditions with the use of sterilized stainless-steel spoons, forceps, spatulas, and scissors.

Containers should be clean, dry, leak-proof, wide-mouthed, sterile, and of a size suitable for samples of the product. Where necessary, containers (metal boxes, cans, bags, packets, etc) should be hermetically sealed. Glass containers should be avoided to minimize the danger of breaking and contaminating the food product. Bags should not be overfilled to avoid puncture by wire closure.

Dry or canned foods that are not perishable and are collected at ambient temperatures need not be refrigerated. Frozen or refrigerated products should be transferred in approved insulated containers to maintain temperatures unchanged. Frozen samples should be collected in pre-chilled containers.

Each sample unit should be identified with a properly marked strip of masking tape or barcode strip. Be aware that the ink of felt pen may penetrate the plastic container.

The person that performed the sampling fills in the sample form in the relevant data fields.

If it is in a paper form, it is sealed in a plastic bag and inserted in the sample packaging. If the form is a spreadsheet it is sent by email to the receiving laboratory.

The time and date of arrival of each sample should be recorded at the laboratory. Frozen samples should be always kept frozen. The containers should be capable of maintaining temperatures of 0-4 °C from collection to arrival at the laboratory for all refrigerated samples, except shellfish and shell stock. Refrigerated products should not be frozen and should be analysed within 36 hours from collection, unless otherwise specified. Shucked shellfish should be kept above freezing but below 10 °C and examined within 6 h of collection but in no case more than 24 h after collection.

As soon as the sample arrives at the laboratory, the analyst should note its general physical condition and record the temperature. Sampling containers should be checked for gross physical defects. Plastic bags and bottles should be inspected for tears, pinholes, and puncture marks. If sample units were collected in plastic bottles, bottles should be checked for fractures and loose lids. If plastic bags were used for sampling, twist wires should not have punctured surrounding bags.

The sample is rejected if:

Sample temperature is not as specified above

Time from sampling is longer than the recommended

Sample container leaks

There are visible cracks or any form of damage to the containers

Sample quantity is less than the recommended

Mandatory fields in the data form are not filled in

If rejected, in the “Reasons to reject the sample” field of the form the exact reason for rejection should be mentioned.

If the sample cannot be analyzed immediately, it should be stored at -20 °C (already frozen samples) or at 0-4 °C (refrigerated unfrozen perishable samples) until examination and for no longer than 36h. Nonperishable, canned, or low-moisture foods should be stored at room temperature until analysis.

Depending on the requested tests, samples are processed according to various ISO methods, or according to recommendations of international or national authorities describing the processing of the sample.

The laboratory analyst completes the "Date and time that laboratory analysis started" field when starting the analysis procedure. If LIMS software is available all metadata are imported.

When results are ready the form is sent by email to the epidemiology department along with the results report, or data may be accessible through LIMS software.

ANNEX I

Official guidelines for pathogen detection and/or enumeration.

For bacterial colony isolation to perform WGS:

ISO 23418:2022 Microbiology of the food chain — Whole genome sequencing for typing and genomic characterization of bacteria — General requirements and guidance

ISO 6579-1:2017/Amd 1:2020 Microbiology of the food chain — Horizontal method for the detection, enumeration and serotyping of *Salmonella* — Part 1: Detection of *Salmonella* spp. — Amendment 1: Broader range of incubation temperatures, amendment to the status of Annex D, and correction of the composition of MSR/V and SC

ISO 6579-1:2017. Microbiology of the food chain — Horizontal method for the detection, enumeration and serotyping of *Salmonella* — Part 1: Detection of *Salmonella* spp.

ISO 6888-1:2021/Amd 1:2023. Microbiology of the food chain — Horizontal method for the enumeration of coagulase-positive staphylococci (*Staphylococcus aureus* and other species) — Part 1: Method using Baird-Parker agar medium — Amendment 1

ISO 6888-2:2021/Amd 1:2023. Microbiology of the food chain — Horizontal method for the enumeration of coagulase-positive staphylococci (*Staphylococcus aureus* and other species) — Part 2: Method using rabbit plasma fibrinogen agar medium — Amendment 1

ISO 7251:2005/Amd 1:2023. Microbiology of food and animal feeding stuffs — Horizontal method for the detection and enumeration of presumptive *Escherichia coli* — Most probable number technique — Amendment 1: Inclusion of performance testing of culture media and reagents

ISO 7932:2004. Microbiology of food and animal feeding stuffs — Horizontal method for the enumeration of presumptive *Bacillus cereus* — Colony-count technique at 30 °C

ISO 10272-2:2017. Microbiology of the food chain — Horizontal method for detection and enumeration of *Campylobacter* spp. — Part 2: Colony-count technique

ISO 10273:2017. Microbiology of the food chain — Horizontal method for the detection of pathogenic *Yersinia enterocolitica*

ISO 11290-1:2017. Microbiology of the food chain — Horizontal method for the detection and enumeration of *Listeria monocytogenes* and of *Listeria* spp. — Part 1: Detection method

ISO 11290-2:2017. Microbiology of the food chain — Horizontal method for the detection and enumeration of *Listeria monocytogenes* and of *Listeria* spp. — Part 2: Enumeration method

ISO 21528-1:2017. Microbiology of the food chain — Horizontal method for the detection and enumeration of Enterobacteriaceae — Part 1: Detection of Enterobacteriaceae

ISO 21528-2:2017. Microbiology of the food chain — Horizontal method for the detection and enumeration of Enterobacteriaceae — Part 2: Colony-count technique

ISO 21567:2004. Microbiology of food and animal feeding stuffs — Horizontal method for the detection of *Shigella* spp.

ISO 15213-2:2023. Microbiology of the food chain — Horizontal method for the detection and enumeration of *Clostridium* spp. — Part 2: Enumeration of *Clostridium perfringens* by colony-count technique

ISO 10272-1:2017/Amd 1:2023. Microbiology of the food chain — Horizontal method for detection and enumeration of *Campylobacter* spp. — Part 1: Detection method — Amendment 1: Inclusion of methods for molecular confirmation and identification of thermotolerant *Campylobacter* spp., the use of growth supplement in Preston broth and changes in the performance testing of culture media

ISO 11866-2:2005. Milk and milk products — Enumeration of presumptive *Escherichia coli* — Part 2: Colony-count technique at 44 °C using membranes

For molecular detection of pathogens by qPCR:

ISO 22174:2005. Microbiology of food and animal feeding stuffs — Polymerase chain reaction (PCR) for the detection of food-borne pathogens — General requirements and definitions

ISO/TS 13136:2012. Microbiology of food and animal feed — Real-time polymerase chain reaction (PCR)-based method for the detection of food-borne pathogens — Horizontal method for the detection of Shiga toxin-producing *Escherichia coli* (STEC) and the determination of O157, O111, O26, O103 and O145 serogroups

ISO/TS 18867:2015. Microbiology of the food chain — Polymerase chain reaction (PCR) for the detection of food-borne pathogens — Detection of pathogenic *Yersinia enterocolitica* and *Yersinia pseudotuberculosis*

ISO 15216-2:2019. Microbiology of the food chain — Horizontal method for determination of hepatitis A virus and norovirus using real-time RT-PCR — Part 2: Method for detection

ISO/TS 17919:2013. Microbiology of the food chain — Polymerase chain reaction (PCR) for the detection of food-borne pathogens — Detection of botulinum type A, B, E and F neurotoxin-producing clostridia

ISO 15216-1:2017. Microbiology of the food chain — Horizontal method for determination of hepatitis A virus and norovirus using real-time RT-PCR — Part 1: Method for quantification

As official molecular detection methods are limited, methods described in ISO documents may be used such as ISO 18744:2016 (Microbiology of the food chain — Detection and enumeration of *Cryptosporidium* and *Giardia* in fresh leafy green vegetables and berry fruits) and instead of microscopy, nucleic acids are isolated from the eluate for molecular detection.

For metagenomic analysis, and/or WGS by selective whole genome amplification, DNA and/or RNA may be isolated directly from the sample.

ANNEX II

Metadata for food samples

	Food sample
1	Sample ID or Barcode*
2	Date and time of sampling*
3	Type of sample/species (e.g. bovine meat, chicken, shells, etc)*
4	Sample condition (frozen, refrigerated, etc)*
5	Additional sample information (e.g., brand name)
6	List of parameters to analyze*
7	Location of sample collection*
8	Link to outbreak (if it is applicable)
9	Date and time of arrival at the laboratory*
10	Sample temperature*
11	Reasons to reject the sample (Yes/No, description)*
12	Other remarks (if any)
13	Laboratory ID*
14	Name of the person that received the sample*
15	Date and time that laboratory analysis started*

Data fields 1-8 are completed by the person performing the sampling, while 9-15 are completed by the laboratory personnel (fields marked with asterisk are mandatory).

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Prikupljanje kliničkih uzoraka opisano je u brojnim dokumentima poput: "J. Michael Miller, Ph.D. Handbook of Specimen Collection and Handling in Microbiology. Second Edition. 1985", "WHO/CDS/CSR/EDC/2000.4. Guidelines for the collection of clinical specimens during field investigation of outbreaks (available at: https://iris.who.int/bitstream/handle/10665/66348/WHO_CDS_CSR_EDC_2000.4.pdf), "Protocol to investigate non-seasonal influenza and other emerging acute respiratory diseases. Geneva: World Health Organization; 2018. (available at: <https://www.who.int/publications/i/item/WHO-WHE-IHM-GIP-2018.2>)" and "Koneman's Color Atlas and Textbook of Diagnostic Microbiology 6th Edition".

Prikupljanje veterinarskih uzoraka opisano je u: "Bildfell R. Collection and submission of laboratory samples from animals. Merck Veterinary Manual" (available at: merckvetmanual.com/clinical-pathology-and-procedures/collection-and-submission-of-laboratory-samples/collection-and-submission-of-laboratory-samples-from-animals), pristupljeno studeni, 2023.

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(WHO Laboratory testing for coronavirus disease (COVID-19) in suspected human cases: interim guidance (<https://www.who.int/publications/i/item/10665-331501>)).

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ISO 9308-1:2014/Amd 1:2016 Water quality - Enumeration of *Escherichia coli* and coliform bacteria - Part 1: Membrane filtration method for waters with low bacterial background flora - Amendment 1

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As official molecular detection methods are limited, methods described in ISO documents may be used such as ISO 18744:2016 (Microbiology of the food chain — Detection and enumeration of

Cryptosporidium and *Giardia* in fresh leafy green vegetables and berry fruits) and instead of microscopy, nucleic acids are isolated from the eluate for molecular detection.

Za metagenomsku analizu, i/ili WGS selektivnom cjelogenomskom amplifikacijom, DNA i/ili RNA se mogu izolirati direktno iz uzorka.

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